

**Delft University of Technology** 

### High-definition microelectrode arrays with scalable, integrated microfluidics in multi-well format for drug screening in a heart-on-a-chip application

Cóndor, Mar ; Braeken, Dries; Rushdy Ahmed, Alaa; Labie, Riet ; Abeille, Fabien ; Meucci, Sandro ; Brandstätter, Birgit ; de Moura, Thiago ; Mastrangeli, Massimo; Stakenborg, Tim

**Publication date** 2023

**Document Version** Final published version

### Citation (APA)

Cóndor, M., Braeken, D., Rushdy Ahmed, A., Labie, R., Abeille, F., Meucci, S., Brandstätter, B., de Moura, T., Mastrangeli, M., & Stakenborg, T. (2023). *High-definition microelectrode arrays with scalable, integrated* microfluidics in multi-well format for drug screening in a heart-on-a-chip application. 263-263. Abstract from 2nd Microphysiological Systems World Summit 2023, Berlin, Germany.

### Important note

To cite this publication, please use the final published version (if applicable). Please check the document version above.

Copyright

Other than for strictly personal use, it is not permitted to download, forward or distribute the text or part of it, without the consent of the author(s) and/or copyright holder(s), unless the work is under an open content license such as Creative Commons.

**Takedown policy** Please contact us and provide details if you believe this document breaches copyrights. We will remove access to the work immediately and investigate your claim.

This work is downloaded from Delft University of Technology For technical reasons the number of authors shown on this cover page is limited to a maximum of 10.



# Abstracts of the 2<sup>nd</sup> Microphysiological Systems World Summit, Berlin, 2023

Volume 11, No. 1 ISSN 2194-0479 doi:10.58847/ap.2301 (2023)

# A LTEAMS

Marcel Leist, Uwe Marx and Peter Loskill **Welcome** 



The 2<sup>nd</sup>

Microphysiological Systems World Summit

26<sup>th</sup>-30<sup>th</sup> June BERLIN 2023 GERMANY Track 1: MPS Development: Bioengineering Models and Readouts

Track 2: MPS for Industrial and Regulatory Application:

Regulatory Application: Standardization, QA, Parallelisation and Automation

Track 3: MPS for Disease Modelling, Safety Testing and Basic Research

Track 4: MPS Highlights Across Disciplines



# - **B**-----

### 496

# High-definition microelectrode arrays with scalable, integrated microfluidics in multi-well format for drug screening in a heart-on-a-chip application

### <u>Mar Cóndor</u><sup>1</sup>, Dries Braeken<sup>1</sup>, Alaa Rushdy Ahmed<sup>1</sup>, Riet Labie<sup>1</sup>, Fabien Abeille<sup>2</sup>, Thiago Moura<sup>3</sup>, Birgit Brandstätter<sup>3</sup>, Massimo Mastrangeli<sup>4</sup>, Sandro Meucci<sup>5</sup> and Tim Stakenborg<sup>1</sup>

<sup>1</sup>imec, Leuven, Belgium; <sup>2</sup>micronit BV, Enschede, The Netherlands; <sup>3</sup>Besi Austria GmbH, Radfeld, Austria; <sup>4</sup>TUDelft, Delft, The Netherlands; <sup>5</sup>Micronit BV, Enschede, The Netherlands

lei.zhang@imec.be

Towards increased throughput and automated workflows for organs-on-a-chip, a novel high-definition electrophysiology multiwell plate is developed in the Moore4Medical project [1]. It consists of an advanced CMOS microelectrode array (MEA) chip with 16 sampling areas, each featuring 1024 electrodes [2], and of a polymeric fluidic cartridge with 16 corresponding microchambers, each connected to an integrated pneumatic peristaltic pump [3]. It can perform 16-well assays with single-cell-resolution electrical recording and integrated microfluidics. Building on an earlier prototype, we present recent developments showing a scalable integration route and excellent fluidic and electrical functionalities of the plate in an easy-to-use system.

Direct bonding of the MEA chip and the cartridge was achieved using accurate glue deposition and a high-precision pick-&-place procedure, resulting in 16 well-aligned, sealed microchambers on top of the MEA. This method is developed for high-volume production. To run the plate, a user-friendly toolbox, which can be placed inside an incubator, is made together with a data acquisition unit, a pneumatic control unit, and custom software. The plate can be clamped directly inside the toolbox, ensuring communication of the MEA with the data acquisition unit by a pogo-pin connection. The pneumatic control unit actuates the 16 integrated pumps simultaneously via three pressure lines. Hence, 16 assays can run in parallel, with simple pipetting steps for cell seeding and media perfusion.

To demonstrate the plate's functions, all the 16 microchambers were primed and filled with water and 1% PBS buffer sequentially. No liquid leakage or bubble entrapment was observed. The voltage scan of the 16 MEA areas showed uniform signals and a clear difference between the two liquids (0.08 mVrms for water; 0.02 mVrms for PBS). Plate biological-validation is currently on-going to study the effect of drugs on electrical and contractile properties of human stem-cell derived cardiomyocytes and to demonstrate the plate potential for OoC standardization and workflow streamlining.

### References

[1] https://moore4medical.eu/organ\_on\_chip

[2] Miccoli, B. et al. (2019). High-density electrical recording and impedance imaging with a multi-modal CMOS multi-elec-

trode array chip. Front Neurosci 13, 641.

[3] Aksoy, B. et al. (2019). Latchable microfluidic valve arrays based on shape memory polymer actuators. *Lab Chip* 19, 608-617.

Presentation: Poster

### 497

## Complementing MPS with mechanistic computer models helps overcome limitations: Translating the drug exenatide from MPS to humans

<u>Oscar Arrestam</u><sup>1,2</sup>, Sophie Rigal<sup>3</sup>, Katharina Schimek<sup>3</sup>, Christian Simonsson<sup>1</sup>, Tilda Herrgårdh<sup>1</sup>, Belen Casas Garcia<sup>1</sup>, Kajsa Kanebratt<sup>4</sup>, Uwe Marx<sup>3</sup>, Liisa Vilén<sup>4,5</sup>, Peter Gennemark<sup>1,4</sup> and Gunnar Cedersund<sup>1,2</sup>

<sup>1</sup>Department of Biomedical Engineering, Linköping University, Linköping, Sweden; <sup>2</sup>SUND sound medical decisions, Linköping, Sweden; <sup>3</sup>TissUse Gmbh, Berlin, Germany; <sup>4</sup>Drug Metabolism and Pharmacokinetics, Research and Early Development, Cardiovascular, Renal and Metabolism (CVRM), BioPharmaceuticals R&D, AstraZeneca, Gothenburg, Sweden; <sup>5</sup>Division of Pharmaceutical Biosciences, Drug Research Program, Faculty of Pharmacy, University of Helsinki, Helsinki, Finland

gunnar.cedersund@gmail.com

Today, realistic organoids can be combined into microphysiological systems (MPS), which are useful for a growing list of applications. Because of this combination, one can study dynamic crosstalk between organs, otherwise only observable in animals. This improved realism creates the potential to replace animal experiments, improve drug development, etc. However, using a purely experimental approach, there are certain limitations that are hard to overcome: the dynamic cross-talk makes data interpretation difficult, and today's MPS still display critical differences to humans, in terms of functionality, volume differences, and missing organs. We therefore propose to combine MPS with mechanistic computer models, which can help overcome these short-comings.

In this presentation, we demonstrate this potential for our two-organ MPS, with liver spheroids and pancreas organ model. Using this MPS, we can study central metabolism, in both healthy conditions and disease states, such as type 2 diabetes and liver steatosis. We develop the computer model by mechanistic hypothesis testing and validate the model by comparing experiments first done in the computer with subsequent MPS results. Using the validated computer model, we can also create more human-like versions of the MPS system. We can, for instance, scale the volumes of the organoids and the circulating media to the human proportions; this makes the consumption of a glucose tolerance test go from the unrealistic 48 h to the typical human timescale of 2-4 h. We can also add already developed computer models for the missing organs – muscle, fat, brain, etc. We demonstrate how this assembling into a