

#### Meta-instrument

### An opto-mechanical platform for imaging near-field optical instruments

Bijster, Roy; Herfst, R.; Klop, W; Hagen, R.; Sadeghian Marnani, Hamed

### **Publication date**

**Document Version** 

Final published version

Citation (APA)

Bijster, R., Herfst, R., Klop, W., Hagen, R., & Sadeghian Marnani, H. (2016). Meta-instrument: An optomechanical platform for imaging near-field optical instruments. Poster session presented at 13th International Workshop on Nanomechanical Sensing, Delft, Netherlands.

Important note

To cite this publication, please use the final published version (if applicable). Please check the document version above.

Copyright

Other than for strictly personal use, it is not permitted to download, forward or distribute the text or part of it, without the consent of the author(s) and/or copyright holder(s), unless the work is under an open content license such as Creative Commons.

Please contact us and provide details if you believe this document breaches copyrights. We will remove access to the work immediately and investigate your claim.

# Meta-Instrument: An Opto-Mechanical Platform for Imaging Near-Field Optical Instruments




Roy Bijster<sup>1,2</sup>, Rodolf Herfst<sup>2</sup>, Wimar Klop<sup>2</sup>, Ronald Hagen<sup>3</sup>, Hamed Sadeghian<sup>2</sup>


- <sup>1</sup>Structural Optimization and Mechanics, Delft University of Technology
- <sup>2</sup>Department of Optomechatronics, Netherlands Organisation for Applied Scientific Research TNO
- <sup>2</sup>Department of Optics, Netherlands Organisation for Applied Scientific Research TNO


Download this poster and the abstract to your phone!

Nanopositioning MEMS device (1) Transparant window, (2) Electronic leads,


(3) Electrodes, (4) Moving plate (SiC), (5) Air gap of 500 nm, (6) Leaf springs









**Fine positioning stage** The stage and interferometers show a large signal coherence up to 10 kHz. This allows for a high bandwidth tracking of the sample.









## Fiber interferometers as distance sensors


Measure the motion of the finepositioning stage. By stabilizing the temperature of the optical path difference, the noise is reduced from 70 nm r.m.s. to 10 nm r.m.s. over 70 kHz.





### **Second generation**

- (1) MEMS nanopositioning device
- (2) Mirror surfaces for interferometers
- (3) Leaf spring that preloads piezo element
- (4) Piezo element for fine positioning
- (5) Fiber interferometer
- (6) Leaf spring that clamps the carbon rod
- (7) Carbon rod of the coarse positioning



## Combined fine and coarse positioning

The new linear actuator can be used for fine positioning and coarse stepping at 0.32 µm/step and velocities up to 3.2 mm/s.

