

## Acoustic Characterization of the CLINicell for Ultrasound Contrast Agent Studies

Beekers, Ines; van Rooij, Tom; van der Steen, Antonius F.W.; de Jong, Nico; Verweij, Martin D.; Kooiman, Klazina

**DOI**

[10.1109/TUFFC.2018.2881724](https://doi.org/10.1109/TUFFC.2018.2881724)

**Publication date**

2018

**Document Version**

Final published version

**Published in**

IEEE Transactions on Ultrasonics, Ferroelectrics, and Frequency Control

**Citation (APA)**

Beekers, I., van Rooij, T., van der Steen, A. F. W., de Jong, N., Verweij, M. D., & Kooiman, K. (2018). Acoustic Characterization of the CLINicell for Ultrasound Contrast Agent Studies. *IEEE Transactions on Ultrasonics, Ferroelectrics, and Frequency Control*, 66(1), 244-246. <https://doi.org/10.1109/TUFFC.2018.2881724>

**Important note**

To cite this publication, please use the final published version (if applicable). Please check the document version above.

**Copyright**

Other than for strictly personal use, it is not permitted to download, forward or distribute the text or part of it, without the consent of the author(s) and/or copyright holder(s), unless the work is under an open content license such as Creative Commons.

**Takedown policy**

Please contact us and provide details if you believe this document breaches copyrights. We will remove access to the work immediately and investigate your claim.

# Acoustic Characterization of the CLINiCell for Ultrasound Contrast Agent Studies

Inés Beekers, *Student Member, IEEE*, Tom van Rooij, *Member, IEEE*,  
 Antonius F. W. van der Steen, *Fellow, IEEE*, Nico de Jong, *Associate Member, IEEE*,  
 Martin D. Verweij, *Member, IEEE*, and Klazina Kooiman, *Senior Member, IEEE*

**Abstract**—Ultrasound contrast agents consist of gas-filled coated microbubbles that oscillate upon ultrasound insonification. Their characteristic oscillatory response provides contrast enhancement for imaging and has the potential to locally enhance drug delivery. Since microbubble response depends on the local acoustic pressure, an ultrasound compatible chamber is needed to study their behavior and the underlying drug delivery pathways. In this study, we determined the amplitude of the acoustic pressure in the CLINiCell, an optically transparent chamber suitable for cell culture. The pressure field was characterized based on microbubble response recorded using the Brandaris 128 ultra-high speed camera and an iterative processing method. The results were compared to a control experiment performed in an OptiCell, which is conventionally used in microbubble studies. Microbubbles in the CLINiCell responded in a controlled manner, comparable to those in the OptiCell. For frequencies from 1 to 4 MHz, the mean pressure amplitude was -5.4 dB with respect to the externally applied field. The predictable ultrasound pressure demonstrates the potential of the CLINiCell as an optical, ultrasound, and cell culture compatible device to study microbubble oscillation behavior and ultrasound-mediated drug delivery.

**Index Terms**—Acoustic characterization, drug delivery, ultra-high speed imaging, ultrasound contrast agents.

## I. INTRODUCTION

ULTRASOUND contrast agents (UCAs) are comprised of gas-filled coated microbubbles (1-10  $\mu\text{m}$  in diameter). Upon ultrasound insonification, the microbubbles compress and expand due to the acoustic pressure. This oscillatory behavior is the characteristic microbubble response used in the clinic for contrast enhancement of diagnostic ultrasound imaging. More recently, their potential to locally enhance vascular drug delivery has been demonstrated [1], [2]. Microbubble oscillation upon ultrasound insonification generates local mechanical stress that can increase endothelial cell membrane permeability by pore formation (i.e. sonoporation), opening of cell-cell junctions, and stimulation of endocytosis [1]. However, the underlying mechanisms of these pathways remain unknown.

Manuscript received June 4, 2018; accepted ..... Date of publication .....; date of current version ..... This work was supported by the Applied and Engineering Sciences TTW (Veni-project 13669), part of NWO.

I. Beekers, T. van Rooij, A. F. W. van der Steen, N. de Jong, M. D. Verweij, and K. Kooiman are with the Department of Biomedical Engineering, Thorax Center, Erasmus MC, Rotterdam, the Netherlands (e-mail: d.beekers@erasmusmc.nl; k.kooiman@erasmusmc.nl).

A. F. W. van der Steen, N. de Jong, and M. D. Verweij are also with the Laboratory of Acoustical Wavefield Imaging, Delft University of Technology, Delft, the Netherlands.

A better understanding of the microbubble's oscillation behavior can aid the design of UCAs for diagnostic ultrasound imaging. In addition, we need to better understand the microbubble-cell interaction to control the different pathways of ultrasound-mediated drug delivery. Microbubble oscillation strongly depends on the ultrasound insonification, composition of shell and gas core, targeting ligands, microbubble size, production method, and their local environment [1]. Therefore, we require a chamber in which microbubble behavior can be characterized experimentally, with or without cellular interaction, which is both ultrasound compatible and optically transparent. If we quantitatively know how the ultrasound pressure field transmits into the chamber, the microbubble can be insonified in a controlled manner. Optical transparency is




Fig. 1. a) Sketch of a CLINiCell with a 25 cm<sup>2</sup> cell culture area. b) Schematic representation (not drawn to scale) of the experimental setup to study microbubble oscillation upon insonification using the Brandaris 128 ultra-high speed camera. The two parallel membranes of the CLINiCell each have a thickness of 50  $\mu\text{m}$ . The microscope was equipped with a 40 $\times$  water immersion objective and a 2 $\times$  lens (Olympus, Tokyo, Japan). The insonification scheme was generated by a waveform generator (Tabor 8026, Tabor Electronics, Tel Hanan, Israel), connected to a 20 dB attenuator (Mini-Circuits, Brooklyn, NY, USA) and a broadband amplifier (ENI A-500, Electronics & Innovation, Rochester, NY, USA).

essential to image microbubble behavior and cellular response with a microscope. Moreover, the chamber should be suitable for cell culture to study ultrasound-mediated drug delivery *in vitro*.

The OptiCell™ (Nunc, Thermo Fisher Scientific, Wiesbaden, Germany) has been extensively used for both microbubble characterization and drug delivery studies [3]. Since the OptiCell is no longer being manufactured, the CLINicell® (MABIO, Tourcoing, France) is an interesting alternative imaging and cell culture chamber. It has recently been used to characterize microbubbles [4], [5] and study ultrasound-mediated drug delivery [6]. However, a thorough acoustic characterization has not been performed yet.

In this study we therefore quantitatively assessed microbubble oscillation in the CLINicell by performing microbubble spectroscopy using the Brandaris 128 ultra-high speed camera. The recordings were iteratively processed using our previously reported method [7].

## II. MATERIAL AND METHODS

The CLINicell (50 μm membrane, 25 cm<sup>2</sup> area) is a cell culture chamber consisting of two parallel gas permeable polycarbonate membranes, with a separation of 5 mm (Fig. 1).

To characterize the acoustic pressure transmitted into the CLINicell *in situ* and non-invasively, microbubbles were used as pressure sensors, as previously described by us [7]. The CLINicell was blocked prior to the experiment with 12 mL of 2% bovine serum albumin (BSA) in phosphate buffered saline (PBS) for 1 h, to prevent unspecific microbubble binding to the membranes. Next, it was rinsed three times with PBS before introducing 12 mL of 10<sup>5</sup> microbubbles/mL. We replicated the experimental set-up, ultrasound insonification parameters, microbubble preparation, and experimental data analysis. Briefly, the DSPC-based microbubbles were insonified at 45° incidence angle with a single-element transducer (1-9 MHz bandwidth, 25 mm focal distance, -6 dB beamwidth at 1 MHz of 1.3 mm, PA275, Precision Acoustics, Dorchester, UK). Microbubble spectroscopy [8] was performed by successively insonifying each individual microbubble while sweeping through a range of transmit frequencies ( $f_T$ ) from 1 to 4 MHz (in steps of 300 kHz). An 8-cycle sine wave burst was applied at 20 kPa external peak negative pressure (PNP), as calibrated in a separate experiment using a 1-mm needle hydrophone. The microbubble oscillation behavior was recorded using the Brandaris 128 ultra-high speed camera [9] (17 million frames/s) combined with a microscope (80x magnification, BX-FM, Olympus, Tokyo, Japan), as depicted in Fig. 1(b).

Using custom-designed image analysis software, microbubble oscillation was quantified as the change in radius over time  $R(t)$ , as previously described [8]. Briefly, the relative excursion was defined as  $x(t) = R(t)/R_0 - 1$ , with resting radius  $R_0$ . The amplitude ( $x_0$ ) of the relative excursion was defined as the maximum after bandpass filtering  $x(t)$  (third-order Butterworth filter, 500 kHz bandwidth, centered at  $f_T$ ). Next, we determined the acoustic pressure amplitude ( $P$ )

experienced by the microbubbles by fitting  $x_0$  in an iterative manner [7] to the harmonic oscillator model

$$x_0 = \frac{|P|/(4\pi^2\rho R_0^2)}{\sqrt{(f_0^2 - f_T^2)^2 + (\delta f_T f_0)^2}} \quad (1)$$

with  $\rho = 10^3$  kg/m<sup>3</sup> being the density of water. The eigenfrequency ( $f_0$ ) of the microbubble is given by

$$f_0 = \frac{1}{2\pi} \sqrt{\frac{1}{\rho R_0^2} \left[ 3\gamma P_0 + \frac{2(3\gamma - 1)\sigma_w}{R_0} + \frac{4\chi}{R_0} \right]} \quad (2)$$

with  $\gamma = 1.07$  the ratio of specific heats for C<sub>4</sub>F<sub>10</sub>,  $P_0 = 10^5$  Pa the ambient pressure,  $\sigma_w = 0.072$  N/m the surface tension, and  $\chi$  the microbubble shell elasticity [8]. The damping coefficient ( $\delta$ ) is given by

$$\delta = \frac{\omega_0 R_0}{c} + 2 \cdot \frac{4\mu}{R_0^2 \rho \omega_0} + \frac{4\kappa_s}{R_0^3 \rho \omega_0} \quad (3)$$

with  $c = 1500$  m/s the speed of sound in water,  $\mu = 10^{-3}$  Pa·s the viscosity of water, and  $\kappa_s$  the microbubble shell viscosity [10]. Therefore, with the iterative scheme not only the pressure experienced by the microbubbles in the CLINicell can be determined, but also the microbubble shell elasticity and viscosity since they are incorporated in  $f_0$  and  $\delta$ . The results were compared to those of control experiments previously performed in an OptiCell, using an identical experimental set-up [10]. All analyses were performed using MATLAB (The MathWorks, Natick, MA, USA).

## III. RESULTS

Single microbubbles with resting radius ( $R_0$ ) from 1.4 to 3.3 μm were recorded upon insonification in the CLINicell (n=38). Fig. 2 shows a typical example of recorded frames and illustrates how the relative excursion amplitude ( $x_0$ ) was determined. This example also demonstrates that bandpass filtering  $x(t)$  removes the low frequency component of compression-only behavior [11], such that the fundamental




Fig. 2. Selected frames of a Brandaris 128 ultra-high speed recording of a single microbubble ( $R_0=2.4$  μm) in the CLINicell, insonified at  $f_T=1.6$  MHz (2 μm scale bar). The dashed line shows the microbubble's relative excursion  $x$  and the black dots indicate the time points of the selected frames. The solid red line shows  $x$  bandpass filtered around  $f_T$  and its maximum, defined as the relative excursion amplitude  $x_0$  (see Fig. 3).




Fig. 3. Relative excursion amplitude ( $x_0$ ) of single microbubbles in the CLINiCell ( $n=38$ , solid red) and OptiCell ( $n=30$ , open blue) as a function of resting radius ( $R_0$ ) when insonified at  $f_T = 1.6$  MHz. The solid data point indicated by the arrow is shown in Fig. 2.




Fig. 4. Normalized pressure amplitude derived from microbubble oscillations in the CLINiCell (solid red) and OptiCell (open/dashed blue). The shaded area corresponds to the interquartile range.

amplitude is obtained. Fig. 3 shows the relative excursion amplitude ( $x_0$ ) as a function of  $R_0$  for every microbubble when insonified at  $f_T = 1.6$  MHz. The microbubble response was similar in both the OptiCell and CLINiCell, with comparable excursion amplitudes and characteristic resonance behavior for  $R_0=2.5$   $\mu\text{m}$  at 1.6 MHz. The pressure experienced by the microbubble ( $P$ ) was obtained by iteratively fitting  $x_0$  to the harmonic oscillator model for 30 iterations. Fig. 4 shows the frequency dependence of  $P$  normalized to the applied external pressure ( $P_A = 20$  kPa PNP). The mean normalized transmitted pressure amplitude from 1 to 4 MHz was similar in both systems:  $-5.4$  dB in the CLINiCell and  $-5.1$  dB in the OptiCell. Moreover, the transmitted pressure showed no clear frequency dependence.

The shell elasticity (with interquartile range between brackets) obtained by fitting the eigenfrequency  $f_0$  (Fig. 5) to equation (2) was similar, with  $\chi = 0.32$  (0.27) N/m in the CLINiCell and  $\chi = 0.32$  (0.25) N/m in the OptiCell. The median shell viscosity obtained from equation (3) was also similar in both systems,  $\kappa_s = 1.1$  (0.5)  $\cdot 10^8$  kg/s in the CLINiCell and  $\kappa_s = 0.7$  (0.4)  $\cdot 10^8$  kg/s in the OptiCell.




Fig. 5. Eigenfrequency ( $f_0$ ) of microbubbles in the CLINiCell (solid red) and OptiCell (open/dashed blue) and the corresponding fit to obtain the shell elasticity. The  $f_0$  of an uncoated microbubble is given by  $\chi = 0$  N/m (black dashed line).

#### IV. DISCUSSION

The acoustic pressure in the CLINiCell was successfully evaluated using microbubbles as non-invasive pressure sensors, confirming the applicability of our iterative processing method [7]. The pressure amplitude transmitted into the CLINiCell was similar to that of the OptiCell and no clear frequency dependence was observed. The microbubbles oscillated in a controlled manner and their shell parameters remained unchanged. We expect an unaltered focal region of the transmitted ultrasound field, since the CLINiCell has two parallel membranes similarly to the OptiCell, in which the focal region remains unchanged [7]. We chose to characterize a CLINiCell with a 50  $\mu\text{m}$  membrane since it is most comparable to the 75  $\mu\text{m}$  membrane of the conventional OptiCell. Also, we assumed that the ultrasound propagation would be affected more by the CLINiCell design with a thicker membrane (125  $\mu\text{m}$ ). The surface area available in the CLINiCell (25  $\text{cm}^2$ ) for cell culture is smaller than in the OptiCell (50  $\text{cm}^2$ ), while the filling volume remains similar (12 mL versus 10 mL). Therefore, microbubble concentration and cell culture protocols might have to be adjusted.

#### V. CONCLUSION

The pressure amplitude inside the CLINiCell was  $-5.4$  dB on average with respect to that of the externally applied pressure, similarly to the OptiCell. Controlled microbubble behavior was achieved in the CLINiCell and the applied pressure can be corrected for the now quantified pressure drop. This demonstrates the potential of the CLINiCell as an alternative chamber to characterize microbubble behavior and to study ultrasound-mediated drug delivery in an optical, ultrasound, and bio-compatible chamber.

#### ACKNOWLEDGMENT

The authors thank Frits Mastik and Robert Beurskens from the Department of Biomedical Engineering, Erasmus MC, for technical assistance during the experiments.

## REFERENCES

- [1] K. Kooiman, H. J. Vos, M. Versluis, and N. de Jong, "Acoustic behavior of microbubbles and implications for drug delivery," *Advanced drug delivery reviews*, vol. 72, pp. 28–48, 2014.
- [2] I. Lentacker, I. De Cock, R. Deckers, S. C. De Smedt, and C. T. W. Moonen, "Understanding ultrasound induced sonoporation: Definitions and underlying mechanisms," *Advanced drug delivery reviews*, vol. 72, pp. 49–64, 2014.
- [3] G. Lajoinie, I. De Cock, C. C. Coussios, I. Lentacker, S. Le Gac, E. Stride, and M. Versluis, "In vitro methods to study bubble-cell interactions: Fundamentals and therapeutic applications," *Biomicrofluidics*, vol. 10, no. 1, 2016.
- [4] H. Shekhar, N. J. Smith, J. L. Raymond, and C. K. Holland, "Effect of Temperature on the Size Distribution, Shell Properties, and Stability of Definity®," 2017.
- [5] A. Nio, A. Faraci, K. Christensen-Jeffries, R. Eckersley, M. Monaghan, J. Raymond, F. Forsberg, and P. Lamata, "The subharmonic amplitude of SonoVue increases with hydrostatic pressure at low incident acoustic pressures," *IEEE International Ultrasonics Symposium, IUS*, no. i, 2017.
- [6] B. H. Lammertink, R. Deckers, M. Derieppe, I. De Cock, I. Lentacker, G. Storm, C. T. Moonen, and C. Bos, "Dynamic Fluorescence Microscopy of Cellular Uptake of Intercalating Model Drugs by Ultrasound-Activated Microbubbles," *Molecular Imaging and Biology*, vol. 19, no. 5, pp. 683–693, 2017.
- [7] I. Beekers, T. van Rooij, M. D. Verweij, M. Versluis, N. de Jong, S. J. Trietsch, and K. Kooiman, "Acoustic Characterization of a Vessel-on-a-Chip Microfluidic System for Ultrasound-Mediated Drug Delivery," *IEEE Transactions on Ultrasonics, Ferroelectrics, and Frequency Control*, vol. 65, no. 4, pp. 570–581, 2018.
- [8] S. M. van der Meer, B. Dollet, M. M. Voormolen, C. T. Chin, A. Bouakaz, N. de Jong, M. Versluis, and D. Lohse, "Microbubble spectroscopy of ultrasound contrast agents," *The Journal of the Acoustical Society of America*, vol. 121, no. 1, pp. 648–656, 2007.
- [9] C. T. Chin, C. Lancée, J. Borsboom, F. Mastik, M. E. Frijlink, N. De Jong, M. Versluis, and D. Lohse, "Brandaris 128: A digital 25 million frames per second camera with 128 highly sensitive frames," *Review of Scientific Instruments*, vol. 74, no. 12, pp. 5026–5034, 2003.
- [10] T. van Rooij, Y. Luan, G. Renaud, A. F. W. van der Steen, M. Versluis, N. de Jong, and K. Kooiman, "Non-linear Response and Viscoelastic Properties of Lipid-Coated Microbubbles: DSPC versus DPPC," *Ultrasound in medicine & biology*, vol. 41, no. 5, pp. 1432–1445, 2015.
- [11] J. Sijl, M. Overvelde, B. Dollet, V. Garbin, N. de Jong, D. Lohse, and M. Versluis, "Compression-only behavior: A second-order nonlinear response of ultrasound contrast agent microbubbles," *The Journal of the Acoustical Society of America*, vol. 129, no. 4, pp. 1729–1739, 2011.