



# FluidBalance

Save Time





Save Money




Increase accuracy



In today's hospitals around 25% of the patients is hospitalized longer than necessary due to complications in their recovery. Many of these complications occur due to wrong treatment plans that are based on inaccurate fluid balance overviews. The FluidBalance system increases accuracy by accurately measuring fluid in- and output on specified times and reduces the workload of the nurse. It automatically updates the fluid balance in the electronic patient files and provides the nurse with reminders to replace or empty a bag.



// The new method for accurate and efficient fluid balance measurements on hospitalized patients. //



The FluidBalance corrects for movement that can influence the measurements before updating the fluid balance data



The FluidBalance hangers are controlled by the main interface to prevent much interaction on non-ergonomic heights




The assistant and the nurse have different levels of access and therefore an own interface to prevent confusion



Drinks and extra losses can be added separately to the FluidBalance system by an assistant or a nurse.

IMPORTANT FEATURES

## CONNECTIVITY SCHEME



Renate Hulst  
FluidBalance  
23-02-2018  
IPD, Medisign specialization

**Committee** S.G van de Geer | R.J.H.G. van Heur (TU Delft)  
D. Vervoorn (Pezy Group)  
R. Pleijhuis | O. Hertgers (TIM Solutions)  
**Company** TIM Solutions

