



SMART DIAGNOSTICS FOR LOW RESOURCE SETTINGS

Target product profiles for devices to diagnose urinary schistosomiasis in Nigeria



Urinary schistosomiasis








The disease can be diagnosed by **counting the** worms' eggs in urine with a microscope.



However, this method requires a **trained operator** and is prone to **human error**.

Smart diagnostic technology



Manual microscopy can be replaced by a device that uses a simple optical system and algorithms to automatically detect Schistosoma eggs in a urine sample. This can offer the following benefits:

- Affordable
- Robust and portable
- Sensitive
- Rapid
- Simple and user-friendly
- Data collection

Diagnostic challenges in Nigeria



There are **no resources** available for diagnosis at primary health care level



There are very **few cases confirmed** in laboratories



Low awareness level amongst community members results in limited care seeking behaviour



Few surveys in the community, so the disease prevalence is unknown



Only children are included in mass drug administration

Promising use cases for smart diagnostics in low resource settings




Mapping of adult populations at risk

Field deployable, affordable diagnostic test, so adult populations at risk can be mapped and treated immediately

Test requirements


- Easy to operate for a lab assistant
- Device is cheaper than a microscope Quantitative and qualitative data is shared with government
- Should comply with WHO guidelines
- Same day result
- Throughput of >50 samples per day
- Should tolerate transport stress
- At least as sensitive and specific as a microscope
- Off grid power supply
- Integrated GPS



Affordable and easy to use diagnostics enable diagnosis at primary healthcare level and sharing of data; referral to laboratories is no longer required

Test requirements

- Easy to operate for a community health worker
- Test should be affordable for community member Confirmed cases should be reported to local
- government • Throughput of 1 sample per day
- Should comply with disease reporting guidelines • At least as sensitive and specific as a microscope
- Should provide treatment dose that has to be administered
- Off grid power supply
- Share data via USB or mobile network



community

Use diagnosis as a tool for health education, to create awareness in the community about the disease

Test requirements

- Easy to operate for a community resource person with very limited test experience
- Visual result, understandable for illiterate
- community members • Test does not make official diagnosis, so positive
- cases should be referred
- Result before the end of sensitization meeting
- Should tolerate transport stress
- Off grid power supply
- Integrated GPS

Merlijn Sluiter Smart Diagnostics for Low Resource Settings 11-05-2020 **Integrated Product Design**

Committee

Dr. ir. J.C. Diehl Prof. dr. ir. J.M.L. van Engelen

